


AQT Series Compact Balances



"The Affordable Choice"

AQT Series Compact Balances

"The Affordable Choice"

The AQT is an economically priced balance tailored for many weighing applications .

- Large and very stable case in modern attractive design
- Stainless steel pan
- Low battery indicator
- Clear backlit LCD display
- Auto power off
- Battery or mains adapter (adapter supplied as standard)
- Various weighing units available depending on model
- Windshield as standard on all models
- Integral security bracket
- Spirit level
- Height adjustable feet

SPECIFICATIONS

Model	AQT 200	AQT 250	AQT 600	AQT 1500	AQT 2000	AQT 5000
Maximum capacity	200g	250g	600g	1500g	2000g	5000g
Readability	0.01g	0.1g		1g		
Resolution	1:20,000	1:2500	1:6000	1:15000	1:2000	1:5000
Repeatability	0.02g	0.1g		0.2g	1g	
Linearity (±)	0.02g	0.1g		0.2g	1g	
Typical stabilisation time	3 Sec	2 Sec		3 Sec	2 Sec	
Units of measure	AQT 200: g,ct,oz,dwt,Taels (Hong Kong, Singapore and Taiwan), AQT 250 and 600: g,ct,oz,dwt, AQT 1500, 2000, 2600 and 5000: kg,g,lb,oz					
Operating temperature	0°C - 40°C					
Power supply	12VAC @ 150mA adapter or 6 x AA batteries					
Calibration	Push button using external mass					
Display	6 digit LCD, 15mm digits with backlight					
Draught shield dimensions	180 x 158 x 82mm / 7.1 x 6.2 x 3.2" (wxdxh)					
Housing	ABS plastic					
Pan size	130mm Ø					
Overall dimensions (wxdxh)	195x240x70mm					
Net weight	1.2kg					
Other features	Low battery, stable and zero indicators					



Optional Accessories:

- In-use Cover



Adam Equipment have offices in the UK, USA, South Africa and China

www.adamequipment.com